

# **Elastomer Testing Services**

The objective of the testing services is to define the basic material properties of elastomeric materials.

# 1. Basic Hyperelastic Properties

3 simple tension specimens, 3 planar tension specimens and 3 equal biaxial specimens are cut from the provided slabs. For foam materials, 3 simple tension specimens, 3 simple shear specimens and 3 simple compression specimens are cut from provided slabs. The specimens are loaded slowly between zero force and a user defined stretch level for 5 loadings and unloadings at up to 4 maximum strain levels, so as to examine the initial stress strain behavior and the "stabilized" stress strain behavior at each of the maximum strain conditions. (5 slabs measuring 150 mm by 150 mm by 1.0 to 2.0 mm thick are required for each temperature condition)

### 2. Volumetric Compression

A specimen is fully constrained and compressed for the purpose of determining the Bulk Modulus of the material. 6.35 mm diameter disks are cut from standard slabs and stacked. The initial slope of this curve is the Bulk Modulus.

#### 3. - 4. Short Term Stress Relaxation Properties

A specimen is stretched to a set strain and held. Stress and time data is collected continuously for a short time. 3 specimens are tested. Simple tension is typically used. This test can be expanded to include multiple strain levels or strain rates.

#### 5. Rate Sensitivity

The rate or speed at which a stress is applied to an elastomer will alter the response of the material. This effect becomes measurable with order of magnitude rate changes.

#### 6. Friction

To measure the aforementioned proportionally factor or coefficient of friction, a 50 mm by 100 mm sled with one material is dragged against a larger second material. Required Materials: a.) Material sheets or plaques greater than 50mm x 150mm for the sled; b.) Materials greater than 100mm x 500mm for the base.

#### 7. Static Tearing

The static tearing experiment is a meaningful way to examine the failure of an elastomer in tearing. A sharp cut is introduced into a planar tension specimen and the specimen is stretched until the cut increases.

## 8. Thermal Properties

Thermal conductivity, thermal diffusivity and specific heat are determined using the transient plane source method.

Thermal expansion is determined by examining the dimensional change in a material specimen TMA.

#### 9. Elastomeric Material Model Calibration Services

The coefficients are determined for an appropriate elastomeric material model and the resulting fit is reported and compared with the experimental test data.

#### 10. Material Preparation Services

Preparing specimens from supplied material slabs is free of charge. To cut specimens from parts or slabs thicker than 2 mm, a fee for skiving services is involved.



# **General Pricing for Elastomer Testing Services**

(Prices are shown in US Dollars)		Room Temp (23C)	-40C < Temp < 200C	37C in Saline
1.	Basic Hyperelastic Properties (typically 5 loadings and unloadings at each of max. 4 strain levels, or single loading)			
	<u>Dense Elastomers</u> 3 Simple Tension Tests 3 Pure Shear Tests 3 Equal Biaxial Tests	TOTAL 1155 210 315 630	TOTAL 1710 315 450 945	TOTAL 2310 420 630 1260
	Foam Elastomers 3 Simple Tension Tests	TOTAL 1155 210	315	n/a
	<ul><li>3 Simple Shear Tests</li><li>3 Simple Compression Tests with Lateral Strain</li></ul>	315 630	450 945	630 420
2.	Volumetric Compression (Bulk Modulus) (3 tests)	250	375	n/a
3.	Standard Short Term Stress Relaxation 3 Simple Tension Tests (2000 seconds, 1 strain level) 3 Fast (10/s) Tension Tests (100 seconds, 1 strain level)	210 440	315 660	420 880
4.	Expanded Short Term Stress Relaxation			
	Dense Elastomers (1Tension, 1 Planar, 1 Biaxial)	1395	2090	2790
	<u>Foam Elastomers</u> (1Tension, 1 Shear, 1 Compression) 5 to 10 strain levels each	1155	1710	2310
5.	Rate Sensitivity Set in Tension (3 tests at 0.01 s <sup>-1</sup> , 5 tests at 1 s <sup>-1</sup> , 5 tests at 50 s <sup>-1</sup> )	2000	3000	n/a
6.	<b>Basic Friction-sled Test</b> (see Friction Test pricing for others) (3 tests at a single pressure less than 0.006 MPa - dry)	230	345	n/a
7.	Static Tearing Experiment (3 tests)	315	450	630
8.	Thermal Properties			
	Thermal Conductivity, Thermal Diffusivity, Specific Heat at 1 temp between –40C and 150C, 3 repetitions250			
	Thermal Conductivity, Thermal Diffusivity, Specific Heat at 5 temps between –40C and 150C, 3 repetitions			600
	Thermal Expansion from –40C to 150C, 3 repetitions_			
9.	9. Elastomeric Material Model Calibration Services			
	Determine Model Coefficients for Hyperelastic Models (i.e. Yeoh, Mooney Rivlin, Ogden)			
10. Material Preparation Services				
	Skiving Parts into Materials price per project			400
ł				

August 28, 2023. Pricing subject to change.

- a. Data is provided in SI units of MPa for stress and non-dimensional strain. The data is delivered via e-mail in an ASCII format.
- b. Customer data and materials will be retained for 1 year after initial data delivery.

Purchase Order, VISA, MasterCard, AMEX, and Discover Card are accepted methods of payment. Terms: NET 30 Days after Delivery of Data